



Notice of Specification Change for EC1260BL, EC1540BL & EC1730BL Combination Warmers:

The maximum temperature set point for the blanket chamber has been modified to 180°F (82°C). All other specifications for this warmer are consistent with the information provided in the Operator's Manual (MN-28646 for 120V units, MN-28672 for 230V units).

Enthermics Medical Systems

ENTHERMICSMEDICALSYSTEMS

An ISO 13485:2016 certified company

PO Box 443, Menomonee Falls WI 53052-0443

W164 N9221 Water St, Menomonee Falls WI 53051

Phone 262-251-8356 | 800-862-9276 | Fax 262-251-7067

generalinfo@enthermics.com

www.enthermics.com

Made in the USA