

WarmWatch®

Temperature recording software



MAINTAIN ACCURATE TEMPERATURE AND PERFORMANCE INFORMATION

ENTHERMICS MEDICAL SYSTEMS

An ISO 13485:2016 certified company

Tel 262-251-8356 | 800-TO-B-WARM

generalinfo@enthermics.com

www.enthermics.com

Safety - ensures temperatures are maintained at proper levels and operating free of any errors

Organization - easily comply with safety standards with fast record retrieval

Time-Saving - automatically records performance parameters



An ISO 13485:2016 certified company

WarmWatch® datalogger included at no charge

ivNext®

- Total capacity: 10 1-liter bags
- Exterior size: 22" H; 18.5" W; 23.5" D
- Datalogger Standard



EC550BL 18" Depth

- Total capacity
 - TriMode chamber: 24 1-liter bottles, 20 1-liter bags or 9-12 blankets
- Exterior size: 25.12" H; 30" W; 21.12" D
- Datalogger Standard



EC550BL 24" Depth

- Total capacity
 - TriMode chamber: 30 1-liter bottles + 36 0.5-liter bottles, 48 1-liter bags or 10-14 blankets
- Exterior size: 25.12" H; 30" W; 28.5" D
- Datalogger Standard



EC1350BL

- Total capacity
 - Blanket Chamber: 24-30 blankets
 - TriMode chamber: 24 1-liter bottles, 20 1-liter bags or 9-12 blankets
- Exterior size: 75" H; 30" W; 21.25" D
- Datalogger Standard

EC1850BL

- Total capacity
 - Blanket Chamber: 30-40 blankets;
 - TriMode chamber: 30 1-liter bottles, 24 1-liter bags or 10-14 blankets
- Exterior size: 75" H; 30" W; 28.5" D
- Datalogger Standard



*Optional feature on DC250L, EC250L & EC350L fluid warmers

**Only use a SanDisk Cruzer Glide USB drive larger than 4 GB