



ENTHERMICS Medical Systems

EC

Mobile Medical Cart

Change your fast-paced patient care traffic patterns with the use of patient carts that roll from one patient to the next, giving you instant access when it matters and decreasing your response time.

Features:

This stylish cart has been engineered to accommodate extra blankets or fluids, or your choice of EC230 blanket warmer or EC230L fluid warmer that easily slides into place in the bottom of the cart. Wrap-up rails make for safe transport of materials carried on top, with rounded corners for easy cleaning. Push handles are integrated for easy maneuvering. The top surface is large enough to accommodate another warmer, allowing you to dispense warmed fluids or blankets as needed. The extended base with full wrap-around vinyl bumper adds stability. Four (4) 5" diameter swivel casters easily maneuver over all floor types — tile, carpet, linoleum, wood and concrete. One (1) caster features a brake for secure positioning.

Specifications:

• **Overall dimensions**

H 41-1/2"
W 23-1/2"
D 30"

• **Top flat surface**

W 20-5/8"
D 24-7/8"

• **Bottom opening**

H 20-1/2"
W 16-1/8"
D 26-1/8"

• **Wrap-up rails double as push handles**

• **Includes one (1) drawer with handle and glides**

H 5-1/2"
W 14-7/8"
D 19-1/2"



Contact your Entermics sales representative for more information.



Corporate Headquarters:

W164 N9221 Water Street • P.O. Box 443
Menomonee Falls, Wisconsin 53052-0443 • U.S.A.



PHONE: 262.251.8356
www.entermics.com

800.TO.B.WARM U.S.A./CANADA
800.862.9276 U.S.A./CANADA

FAX: 262.251.7067
800.329.8744